

Zeiss Stemi 508

mit PL 10x 23 Br. foc. Okular | *with PL 10x 23 Br. foc. eyepiece*

Vorsatzoptik <i>Front optic</i>	Arbeitsabstand <i>Free working distance</i>		Objektfeld <i>Object field</i>
0,3x	287 mm	1,9x – 15x	122,7 – 15,3 mm
0,4x	211 mm	2,5x – 20x	92 – 11,5 mm
0,3-0,5x	230 – 90 mm	1,9x – 25x	122,7 – 9,2 mm
0,63x	127 mm	3,9x – 31,5x	58,4 – 7,3 mm
1x ohne Vorsatzoptik <i>without front optic</i>	92 mm	6,3x – 50x	36,8 – 4,6 mm
1,5x	53 mm	9,4x – 75x	24,5 – 3,1 mm
2x	35 mm	12,5x – 100x	18,4 – 2,3 mm

Zeiss Stemi 508

mit PL 16x 14 Br. foc. Okular | *with PL 16x 14 Br. foc. eyepiece*

Vorsatzoptik <i>Front optic</i>	Arbeitsabstand <i>Free working distance</i>		Objektfeld <i>Object field</i>
0,3x	287 mm	3x – 24x	85,3 – 10,7
0,4x	211 mm	4x – 32x	64 – 8
0,3-0,5x	230 – 90 mm	3x – 40x	85,3 – 6,4
0,63x	127 mm	6,3x – 50x	40,6 – 5,1
1x ohne Vorsatzoptik <i>without front optic</i>	92 mm	10x – 80x	25,6 – 3,2
1,5x	53 mm	15x – 120x	17,1 – 2,1
2x	35 mm	20x – 160x	12,8 – 1,6